





Data
Science


PVC Detection

Premature Ventricular Contractions detection

Projecto de Engenharia e Ciencia de Dados



Premature Ventricular Contraction (PVC)



Summary


■ 1 | Objectives

- ECG - Electrocardiogram
- PVC detection
- Datasets / validation
- Requirements

1. Objectives

■ Clinical problem

- Context
 - **PVC** – premature ventricular contractions
 - A specific type of arrhythmia
- Main goal
 - Given an electrocardiogram, **how many PVC/ hour** ?
 - In other words, classify an **ECG beat as: normal or abnormal {0,1}**



Summary

- Objectives
- **2| ECG - electrocardiogram**
- PVC detection
- Datasets / validation
- Requirements

2. Eletrocardiogram



The heart's function is to provide oxygen to the human body organs .

How it works ?

Myocardium, valves, atria, ventricles, ...

Large circulation / small circulation

...



2. Eletrocardiogram

■ Electrocardiogram


- Measure the electrical activity of the heart





2. Eletrocardiogram

■ Cardiac cycle – electrical activity

- Sinoatrial node (SA)
- Atrioventricular node (AV)
- Purkinje fibers




- Systolic = contraction
- Diastolic = relaxation



2. Eletrocardiogram

■ Cardiac cycle


- P and T waves
- QRS complex
- R peak



2. Eelectrocardiogram

■ Normal rhythm


- 60-100 beats / min
- P waves before the QRS complex
- PR interval [0.12 - 0.20] seconds, and approximately constant
- QRS interval [0.06 - 0.12]



2. Eelectrocardiogram

■ Arrhythmias

- Types of arrhythmias ??
- Several types (dozens)!
- **Two main types**
 - Rhythm (regular ECG)
 - Morphology (irregular ECG)




Summary

- Objectives
- R peak detection
- **3| PVC detection**
- Datasets / validation
- Requirements

3. PVC detection

■ Sinoatrial Node

- A natural pacemaker




© MAYO FOUNDATION FOR MEDICAL EDUCATION AND RESEARCH. ALL RIGHTS RESERVED.

3. PVC detection

■ Premature ventricular contractions

- Ectopic focus !



3. PVC detection

■ PVC – Premature ventricular contraction

- A **focus in the ventricles** can "trigger" and originate an extra ventricular systole.




- A single PVC is not clinically significant
- However, PVC may progress to dangerous conditions (ventricular fibrillation)

3. PVC detection

■ How to detect a PVC

- The beat occurs before the expected
- The complex QRS presents an excessive width/duration, when compared with a normal beat




Summary

- Objectives
- R peak detection
- PVC detection
- **4 | Datasets / validation**
- Requirements

4. Data sets / validation

■ DATPVC (360 Hz) - PVC

- DAT.ecg - signal
- DAT.ind - R peaks indexes
- DAT.pvc - 0/1 - normal / PVC



- 11 records / patients
- Approximately 30 minutes each
- Sampling rate 360 Hz


Summary

- Objectives
- R peak detection
- PVC detection
- Datasets / validation
- **5|** Requirements

5. Requirements

■ For each patient (file)

- Given an ECG register



- **Compute / visualize**

- The number of total beats
- The mean heart rate
- The number of PVCs
- The number of PVC/hour
- The F1 score of the PVCs classification system
- The average ECG (**figure**) using normal cycles

Example

> 220 beats

> 65 bpm

> 39

> 14

> 0,78

